



Intra-sample three-point quantitative calibration

using multiple stable isotope labelled QconCATs

Florian Christoph Sigloch, PolyQuant GmbH, Bad Abbach

www.polyquant.com

Introduction

Quantitative concatamers (QconCATs) have been used since almost 20 years as quantification standards for multi-protein panels in targeted and untargeted proteomics assays. Classical approaches use heavy stable isotope labelled QconCATs for single-point internal calibration. The most common labelling strategy uses labelled arginine and lysine residues, a technique known as SILAC, which is also often used for peptide standards.


Peptides are quantified by adding the labelled standards in known amounts to the samples and by comparing the signal intensity of unlabelled, native peptides with their labelled counterparts.

Here, we show the advantages and drawbacks of metabolic QconCAT labelling using ¹⁵N salt and/or ¹³C glucose compared to classical SILAC labelling. By spiking multiple differently labelled QconCATs at increasing concentrations into the same sample, we establish an intra-sample three-point calibration curve.

Concept

QconCAT technology


QconCATs are artificial concatamers of quantitative peptides for several target proteins. They are typically used as internal standards in LC-MS/MS based proteomics, but can also serve as external standards. QconCATs can be produced as recombinant proteins in *E.coli*. For heavy isotope labelling, bacterial cultures are provided with either heavy labelled amino acids (SILAC), ¹⁵N salt and/or ¹³C glucose.



Three-point internal calibration

MS signal intensities show a linear positive correlation with protein / peptide abundance. However, the slope of the line often differs from 1, i.e. a doubling of the peptide's intensity does not lead to a doubling of the observed signal intensity. This phenomenon reduces the accuracy of one-point calibrations, especially when the signal intensities of standard and unknown differ strongly.

Therefore, quantitative clinical MS assays must often use external calibration curves. A drawback of external calibration is lacking intra-sample controls that can correct e.g. for matrix effects or ion suppression. Also, external calibration curves need additional measurements and thus increase the needed machine time.



Results

Analysis of labelling strategies

A common problem in quantitative proteomics is incomplete labelling of the standards during production. We discriminate two sources of error that may disturb exact quantification:

1) Contamination with unlabeled peptides

Standards are often contaminated by fully unlabelled peptides that have the same mass than native proteins. The resulting signals in the light trace cannot be distinguished from the native protein and thus limit the linearity and sensitivity of the method. Unlabelled peptides occur in relevant amounts only in SILAC and ¹³C labelling, the signal detected in ¹⁵N and ¹³C¹⁵N labelling is negligible.




Figure 1: Percentage of unlabeled peptide signal of total peptide signal. SILAC peptides show a median unlabeled signal of 1.54 %, ¹³C labelled peptides a median of 0.16 %. The unlabeled signal in ¹⁵N and ¹³C¹⁵N peptides is less than 0.01% and can be mostly attributed to noise.

2) Partial labelling

A different problem arises when not all ¹²C or ¹⁴N atoms of a given peptide are replaced by heavy isotopes. Incomplete labelling results in the observation of peaks that are lighter than the expected monoisotopic peak. While ¹⁵N and SILAC labelling show an isotope pattern comparable to unlabelled peptides, ¹³C and ¹³C¹⁵N labelling show a shift towards lighter isotopes.




Figure 3: Precursor isotopes observed in ¹³C labelling. The percentage of not fully labelled standard increases with the peptide length, while the percentage of the monoisotopic precursor drops.




Figure 2: Isotope distribution. Isotopic peaks were identified and quantified using Skyline. All peaks with at least 1% signal compared to the monoisotopic precursor were included.

Skyline XICs





Figure 4: Extracted ion chromatograms from labelled QconCATs. The same QconCAT was expressed five times with different labelling: unlabelled, SILAC-¹³C¹⁵N @ Lys & Arg, metabolic ¹³C, metabolic ¹³C¹⁵N and metabolic ¹⁵N. Each QconCAT was measured separately and the XICs for the light and heavy peptide LISEVIGER was extracted using the Skyline software.

Proof of principle

Differentially labelled QconCATs were spiked into the same sample in increasing amounts, to calculate an intra-sample calibration curve. QconCATs with ¹³C, ¹³C¹⁵N and ¹⁵N were spiked in a 2 : 10 : 50 ratio. Unlabelled QconCAT served as an “unknown” and was spiked at a ratio of 5. All signals were corrected for the observed differences in the isotope distribution and for slight differences in the peak integration.



Label	Signal	True conc	Calculated
13C	6.6x 10 ⁷	2	-
13C15N	3.5x 10 ⁸	10	-
15N	1.8x 10 ⁹	50	-
unlabelled	1.7x 10 ⁸	5	4.997

Figure 5: Intra-sample linear regression. A linear fit was calculated from the different heavy signals. The unlabelled QconCAT was used as an unknown. Its concentration was back-calculated based on its signal from the linear regression.

Conflict of interest

Florian Christoph Sigloch is employed by PolyQuant GmbH

